

Automation of seizure detection using EEG signals.

Prepared by:

Abir LAKHAL

Mehdi HMAIDI

Racha MAGHERBI

Khaled CHARAABI

Nada BELAIDI

Fares KHECHANA

Presentation Plan :

01.

Introduction

02.

Business & Data
Comprehension

03.

Data Preparation

04.

Data Modeling &
Evaluation

05.

Deployment

06.

Product Advertising

07.

Conclusion



01

Introduction



01.1

Problem Identification

Problem :

- **Epilepsy is a neurological disease that results in abnormal electrical activity in the brain.**
- **A sudden increase in electrical activity in the brain**
- **A temporary disruption of communication between neurons.**



Detection :

Analysis of EEG (electroencephalogram) signals to detect and categorize disease trends:

- **Long and laborious visual examination.**
- **Lots of hours to review the recording.**
- **A heavy burden on neurologists and reduces their efficiency.**



Target audience :

Neurologists

&

Neurology Institutes.

Target audience' needs :

**An application that helps
neurologists detect epileptic
seizures is needed.**






01.2

Product Description

General Description :

Application dedicated to neurologists and neurology departments in hospitals.

The application is based on the detection, clinically, of epileptic seizures, i.e. scalp electroencephalogram (EEG) analysis.



Product-Service Description :

Help the neurologist to :

- Properly detect and categorize the disease of epilepsy.**
- Keep the efficiency of neurologists.**
- Avoid medical errors.**
- Increase the diagnosis accuracy.**
- Cost reduction : decrease the number of experts in hospitals.**
- Save time: reduce the need of an expert.**
- Reduce the EEG training time**


Promises and potential :

**Classification of EEG signals:
epileptic and non-epileptic.**

**Processing of MRI
images to detect the
epileptic area.**

**Processing of medical
reports to classify
epileptic of epilepsy.**

**Voice recognition of
doctors and automation
of writing medical reports
(voice to text).**






02.

Business & Data Comprehension

Business Understanding :



Description of the specific objective selected

Create an epilepsy detection platform dedicated to neurologists in order to help them detect and categorize the disease in a short time

Data Science Objectives :

- **Creating an epilepsy detective system.**
 - **Facilitating neurologists' jobs through our system.**
 - **Detecting much more cases.**
 - **Modernize the Tunisian medical field.**
- **Identifying the suitable technologies for our business objectives.**
- **Training and deploying fast and efficient Deep Learning models.**

Key Results :

- **Using NLP for text mining.**
- **Using LSTM for EEG classification.**
- **Using machine learning algorithms for text classification.**
- **Using R-CNN for MRI analysis.**
- **Implementing the API in a demo web application (Flask app).**

Data Science Project Steps :

Feature engineering

Data understanding.
Data cleaning.
Creating new features for the database.

External data

- **Add an EDF file that describes the brain signal of the patient.**
- **Use an MRI Dataset.**
- **Use medical reports to classify epilepsy type.**

Data modeling

Classification.
Scoring.


Deployment

Creating a web platform for this epilepsy detective system (Using Flask).



03.

Data Preparation



03.1

Internal Data Preparation

EEG SIGNALS :

THE EEG signals were recorded using a standard 10 to 20 channels. The complete data consists of five sets (A to E), each containing 100 instances .

Each file is a recording of brain activity for 23.6 seconds. The time series corresponding is sampled at 4097 data points. Each data point is the value of the EEG recording at a different time. So we have a total of 500 individuals each with 4097 data points.

- 1- Loading files from folders then we extracting the data into 5 small tables, each for every set, then we assembled all of data into a big dataframe.**
- 2- Checking for missing data.**
- 3- Changing the labels to binary (0 for healthy, 1 for epileptic).**
- 4- Normalizing the data.**
- 5- CHI2 statistical test (for feature selection).**

```
###Statistical Test to determine whether input features are relevant to the outcome to be predicted.
```



```
###P-value <= 0.05 significant result
```

```
###P-value > 0.05 not significant result
```

Out[45]:

	file_name	1	2	3	4	5	6	7	8	9	...	4089	4090	4091	4092	4093	4094	4095	4096	4097	label
0	Z001.txt	12.0	22.0	35.0	45.0	69.0	74.0	79.0	78.0	66.0	...	-28.0	-21.0	-14.0	-14.0	-25.0	-28.0	-11.0	8.0	77.0	A
1	Z002.txt	-56.0	-50.0	-64.0	-91.0	-135.0	-140.0	-134.0	-114.0	-115.0	...	-82.0	-114.0	-138.0	-159.0	-172.0	-180.0	-173.0	-162.0	-82.0	A
2	Z003.txt	-37.0	-22.0	-17.0	-24.0	-31.0	-20.0	-5.0	14.0	31.0	...	-52.0	-23.0	-14.0	-5.0	-3.0	7.0	3.0	4.0	82.0	A
3	Z004.txt	-31.0	-43.0	-39.0	-39.0	-9.0	-5.0	18.0	7.0	-12.0	...	-32.0	-40.0	-23.0	-1.0	11.0	12.0	-6.0	10.0	33.0	A
4	Z005.txt	14.0	26.0	32.0	25.0	16.0	8.0	8.0	12.0	11.0	...	-19.0	-29.0	-35.0	-51.0	-55.0	-58.0	-32.0	-6.0	-17.0	A
...
95	S096.txt	-40.0	-58.0	-75.0	-88.0	-89.0	-81.0	-67.0	-52.0	-28.0	...	32.0	32.0	18.0	6.0	-3.0	-10.0	-13.0	-16.0	-151.0	E
96	S097.txt	187.0	44.0	-147.0	-368.0	-550.0	-657.0	-665.0	-581.0	-442.0	...	510.0	562.0	607.0	667.0	748.0	763.0	703.0	446.0	-537.0	E
97	S098.txt	-438.0	-561.0	-622.0	-581.0	-460.0	-295.0	-164.0	-70.0	3.0	...	443.0	399.0	319.0	196.0	40.0	-47.0	-118.0	-163.0	-56.0	E
98	S099.txt	-476.0	-518.0	-521.0	-362.0	-68.0	175.0	289.0	184.0	15.0	...	-261.0	-248.0	-147.0	36.0	224.0	299.0	246.0	556.0	276.0	E
99	S100.txt	23.0	144.0	228.0	260.0	255.0	218.0	178.0	126.0	60.0	...	-127.0	-123.0	-152.0	-231.0	-272.0	-272.0	-155.0	6.0	-221.0	E

500 rows x 4099 columns



Clinical reports :

The dataset contains 522 clinical report texts, containing information about each patient and his treatment, divided in epileptic and none epileptic.

Data Preparation Steps :

- Collecting of clinical text reports.**
- Appending of all readlines in one matrix.**
- Applying the natural language processing.**
- Converting text to lowercase.**
- Removing punctuation.**
- Removing whitespaces.**
- Sentence segmentation.**
- Stopwords removal.**
- Count vectorizing.**
- word level tf-idf**






03.2

External Data Preparation

MRI :

The dataset contains 50 brain T1-weighted MRI volumes with hippocampus labels (13 males / 27 females).

Age range : 15-64.




We have divided the dataset into two groups of training and testing (50% each). We used the libraries nilearn and nibabel to visualize and understand our data. Our data was composed of .HDR files so we had to extract a 3D matrix from each file and flattened them into one array and then normalized it.

The preparation of our data was based on :

1-Putting all our MRI files into the same shape :


Since our dataset files have different shapes,sizes and dimensions , we will be reducing all of them to one main size.

2- Smoothing them to become more clear



3- Cropping them to get rid of all the unnecessary voxels (voxels outside the brain) using masks .

4- Creating a mask for the hippocampus MRI labels and then cropping them too.



VOICE RECOGNITION API :


**Speech recognition helps us save time by speaking instead of typing.
We used google speech to text API to help medical stuff and doctors record
their clinical reports without needing a pen.**





04.

Data Modeling & Evaluation



04.1

Internal Data Modeling & Evaluation

EEG SIGNALS :

To find out which is the most suitable model for our internal data we've started by applying SIMPLE RNN model.
We came out with 0.75 as a training accuracy for this model.

```
13/13 [=====] - 0s 3ms/step - loss: 0.6755 - val_loss: 0.7820
Epoch 93/100
13/13 [=====] - 0s 4ms/step - loss: 0.7126 - val_loss: 0.7821
Epoch 94/100
13/13 [=====] - 0s 4ms/step - loss: 0.7266 - val_loss: 0.8117
Epoch 95/100
13/13 [=====] - 0s 4ms/step - loss: 0.7459 - val_loss: 0.8153
Epoch 96/100
13/13 [=====] - 0s 3ms/step - loss: 0.6796 - val_loss: 0.8428
Epoch 97/100
13/13 [=====] - 0s 3ms/step - loss: 0.6867 - val_loss: 0.8323
Epoch 98/100
13/13 [=====] - 0s 3ms/step - loss: 0.7201 - val_loss: 0.7561
Epoch 99/100
13/13 [=====] - 0s 3ms/step - loss: 0.7195 - val_loss: 0.7794
Epoch 100/100
13/13 [=====] - 0s 3ms/step - loss: 0.7137 - val_loss: 0.7542
4/4 [=====] - 0s 997us/step - loss: 0.7542
```

Out[19]: 0.7542365193367004


Then we've chosen to apply LSTM model that provided us with a better training accuracy (0.95).

```
Entrée [24]: from sklearn.metrics import accuracy_score
pred = model.predict(x_test)
predict_classes = np.argmax(pred,axis=1)
expected_classes = np.argmax(y_test,axis=1)
print(expected_classes.shape)
print(predict_classes.shape)
correct = accuracy_score(expected_classes,predict_classes)
print(f"Training Accuracy: {correct}")

(100,)
(100,)
Training Accuracy: 0.95
```

Finally we tried some basic Machine learning models.

As a result, 3 models (RandomForest, ExtraTrees, GradientBoosting) have shown an overfitting and the rest (SVM, DecisionTree, KNN, LogisticRegression) came out with different accuracy with the best one is 0.85 for Decision Tree .



Clinical Reports :

```
[87] #naive_bayes
gnb = naive_bayes.MultinomialNB()
gnb.fit(x_train,y_train)
predictions = gnb.predict(x_test)
metrics.accuracy_score(predictions, y_test)
```

0.8952380952380953

```
▶ knn = KNeighborsClassifier(2)
knn_model = knn.fit(x_train, y_train)
y_pred_knn =knn_model.predict(x_test)

print('Accuracy of K-NN classifier on training set: {:.2f}'
      .format(knn.score(x_train, y_train)))
print('Accuracy of K-NN classifier on test set: {:.2f}'
      .format(knn.score(x_test, y_test)))
```


```
↳ Accuracy of K-NN classifier on training set: 0.96
Accuracy of K-NN classifier on test set: 0.91
```



Here we used other machine learning models (CART, Linear model, Random Forest) but the result was an overfitting. For the rest we had a better accuracy with the KNN model.

```

CrossValMeans  CrossValerrors  Algorithm
0      1.000000      0.000000  Random forest
1      1.000000      0.000000  linear model
2      0.825761      0.057021  naive bayes
3      1.000000      0.000000  XGBC
4      0.838549      0.048592  KNN
5      0.991414      0.010176  CART
/usr/local/lib/python3.7/dist-packages/seaborn/_decorators.py:43: FutureWarning: Pass t
FutureWarning

```





04.2

External Data Modeling & Evaluation

MRI :

Since 3D U-Net model is the most suitable for our Data ,we've chosen to apply it on our MRI files.

We proceeded to use MONAI and Catalyst, two PyTorch-based frameworks in order to implement our model since it wasn't available anywhere and had to be coded from scratch.

```
# create UNet, DiceLoss and Adam optimizer  
# device = torch.device("cuda:0") # you don't need device, because Catalyst uses autoscaling  
model = UNet(  
    dimensions=3,  
    in_channels=1,  
    out_channels=1,  
    channels=(16, 32, 64, 128, 256),  
    strides=(2, 2, 2, 2),  
    num_res_units=2,  
)
```

This was our training output, we had 0.94 as an accuracy value.

```
ice_metric=0.9410 | loss=-7.309e-01
2/2 * Epoch 2 (valid): _timer/_fps=0.2591 | _timer/batch_time=6.4243 | _timer/data_time=0.2666 | _timer/model_time=6.1577 | d
ice_metric=0.9441 | loss=-7.498e-01
[2021-04-29 00:54:11,730]
2/2 * Epoch 2 (_base): lr=0.0010 | momentum=0.9000
2/2 * Epoch 2 (train): _timer/_fps=0.1138 | _timer/batch_time=9.6863 | _timer/data_time=6.5672 | _timer/model_time=3.1191 | d
ice_metric=0.9410 | loss=-7.309e-01
2/2 * Epoch 2 (valid): _timer/_fps=0.2591 | _timer/batch_time=6.4243 | _timer/data_time=0.2666 | _timer/model_time=6.1577 | d
ice_metric=0.9441 | loss=-7.498e-01
INFO:metrics_logger:
2/2 * Epoch 2 (_base): lr=0.0010 | momentum=0.9000
2/2 * Epoch 2 (train): _timer/_fps=0.1138 | _timer/batch_time=9.6863 | _timer/data_time=6.5672 | _timer/model_time=3.1191 | d
ice_metric=0.9410 | loss=-7.309e-01
2/2 * Epoch 2 (valid): _timer/_fps=0.2591 | _timer/batch_time=6.4243 | _timer/data_time=0.2666 | _timer/model_time=6.1577 | d
ice_metric=0.9441 | loss=-7.498e-01
Top best models:
C:\Users\msi\AppData\Local\Temp\tmpgq7413gq\logs\checkpoints/train.2.pth          0.9441
=> Loading checkpoint C:\Users\msi\AppData\Local\Temp\tmpgq7413gq\logs\checkpoints\best_full.pth
loaded state checkpoint C:\Users\msi\AppData\Local\Temp\tmpgq7413gq\logs\checkpoints\best_full.pth (global epoch 2, epoch 2,
stage train)
```

VOICE RECOGNITION API :

```
import speech_recognition as sr
```

```
r = sr.Recognizer()
```

```
harvard = sr.AudioFile('enree.wav')  
with harvard as source:  
    r.adjust_for_ambient_noise(source)  
    audio = r.record(source)
```

```
r.recognize_google(audio)
```

```
'today I want you to ask yourself this one question why not you why not you to do something for work that you love'
```



05.

Deployment

**We used Flask, which is an open-source micro framework
for web development in Python.**



Flask

web development,
one drop at a time



EPILEPTICA

Welcome to our website

invisible seizures now visible



06.

Product Advertising

EPILEPTICA:

- An application that helps neurologists detect epileptic seizures.
- Increases the diagnosis accuracy.
- Avoids medical errors.
- Cost reduction.
- Saving time.

4P: Marketing Mix



- Online.
- Google play.
- App store.


- Advertising.
- Public relations.
- Personal sales.
- Sales promotion.
- Online Marketing.
- Social media.
- Sponsorship.

- Perception of values.

Logo :



EPILEPTICA




EPILEPTICA



EPILEPTICA




EPILEPTICA




EPILEPTICA


Graphic Chart :



ec7c90



df206a



ec7c90




The background features decorative circuit board patterns in the corners, rendered in a light gray color. These patterns consist of interconnected lines and small circles, resembling a stylized electronic circuit.

Marketing materials.

Business Card :



Poster :



Call-to-action :





07.

Conclusion

Our application will facilitate neurologists' job face to epilepsy disease, in terms of time saving and cost reduction.



EPILEPTICA

Contact us :



www.facebook.com/epilepticaApp



[@epilepticaApp](https://www.instagram.com/epilepticaApp)



epilepticaapp@gmail.com



**Thanks
for your
attention.**